



COUNTRY ALU

1416SR & 1417SR

SILLA DE ALUMINIO PLEGABLE
ALUMINUM FOLDING WHEELCHAIR



LIGERA
LIGHTWEIGHT



HASTA 125KG.
UP TO 125KG.



ALUMINIO
ALUMINIUM

1416SR

1417SR

AZUL OSCURO
DARK BLUE

MEDIDAS
MEASURE

	41		55										
1416SR	43	100	58									12	
	45		61										
	47		64	30	96	51	43	42	125				
1417SR	41	110	57								15		
	43		60										
	45		63										
	47		66										
	cm.											Kg.	

Estructura Frame	Aluminio, plegable. Aluminum, folding.	
Ruedas Delanteras Front Wheels	180mm macizas. 180mm PU.	
Ruedas Traseras Rear Wheels	315mm macizas. 315mm PU.	600mm macizas. Extracción rápida. 600mm PU. Quick Release
Reposabrazos Armrests	Abatible y extraíble. Folding and detachable.	
Reposapiés Footrests	Abatible interior-exterior. Altura regulable. Swing in and out footrest. Height adjustment	
Respaldo Backrest	No plegable. Bolsillo portaobjetos. Non-foldable. Back pocket.	
Accesorios Accessories	Anti-vuelcos. Reposapiés elevables. Porta oxígeno. Porta suero. Reposacabezas. Porta bastones. Anti-tipping wheels. Elevating footrests. Oxygen holder. Serum holder. Headrest. Cane holder.	



AZUL CLARO
LIGHT BLUE

MEDIDAS
MEASURE

	41		59,5										
1418SR	43	98	61,5	29								14 - 15	
	45		63,5										
	47		65,5		92	48,5	43	42	125				
1419SR	41	90	58	27							12 - 13		
	43		60										
	45		62										
	47		64										
	cm.											Kg.	

FORUM

1418SR & 1419SR

SILLA DE ALUMINIO PLEGABLE DOBLE CRUCETA
DOUBLE CROSS BAR ALUMINUM FOLDING WHEELCHAIR



LIGERA
LIGHTWEIGHT



HASTA 125KG.
UP TO 125KG.



ALUMINIO
ALUMINIUM



DOBLE CRUCETA
DOUBLE CROSSBAR

1418SR

1419SR

Estructura Frame	Aluminio, plegable. Doble cruceta. Aluminum, folding. Double crossbar.	
Ruedas Delanteras Front Wheels	200mm macizas. 200mm PU.	
Ruedas Traseras Rear Wheels	300mm macizas. 300mm PU.	600mm macizas. Extracción rápida. 600mm PU. Quick Release
Reposabrazos Armrests	Abatible y extraíble. Folding and detachable.	
Reposapiés Footrests	Abatible interior-exterior. Altura regulable. Swing in and out. Height adjustable.	
Respaldo Backrest	No plegable. Non-foldable.	
Accesorios Accessories	Anti-vuelcos. Reposapiés elevables. Porta oxígeno. Porta suero. Reposacabezas. Porta bastones. Anti-tipping wheels. Elevating footrests. Oxygen holder. Serum holder. Headrest. Cane holder.	