



RAMBLA

1404SR & 1405SR

SILLA DE ACERO PLEGABLE
FOLDING STEEL WHEELCHAIR



HASTA 125KG.
UP TO 125KG.



ACERO
STEEL

1405SR

1404SR

PLATEADO
SILVER

MEDIDAS
MEASURE



1404SR	41	101	59	90	43	49	43	125	18,5
	43		61						
	45		63						
1405SR	41	105	56	30	96	52	43	41	16
	43		58						
	45		60						
cm.									Kg.

Estructura Frame	Acero, plegable. Steel, folding.	
Ruedas Delanteras Front Wheels	200mm macizas. 200mm PU.	
Ruedas Traseras Rear Wheels	300mm macizas. 300mm PU.	600mm macizas. Extracción rápida. 600mm PU. Quick Release
Reposabrazos Armrests	Abatible y extraíble. Folding and detachable.	
Reposapiés Footrests	Abatible exterior. Altura regulable. Swing away. Height adjustable.	
Respaldo Backrest	No plegable. Bolsillo portaobjetos. Non-foldable. Back pocket.	
Accesorios Accessories	Anti-vuelcos. Reposapiés elevables. Porta oxígeno. Porta suero. Reposacabezas. Porta bastones. Anti-tipping wheels. Elevating footrests. Oxygen holder. Serum holder. Headrest. Cane holder.	



COUNTRY

1412SR & 1413SR

SILLA DE ACERO PLEGABLE
FOLDING STEEL WHEELCHAIR



HASTA 135KG.
UP TO 135KG.



ACERO
STEEL

1412SR

1413SR

AZUL OSCURO
DARK BLUE

MEDIDAS
MEASURE



1412SR	39,5	99	53	26	96	52	43	41	135	17
	42		55							
	44		57							
	46		59							
1413SR	39,5	105	58	30	96	52	43	41	135	18
	42		60							
	44		62							
	46		64							
cm.									Kg.	

Estructura Frame	Acero, plegable. Steel, folding.	
Ruedas Delanteras Front Wheels	180mm macizas. 180mm PU.	
Ruedas Traseras Rear Wheels	315mm macizas. 315mm PU.	600mm macizas. Extracción rápida. 600mm PU. Quick Release
Reposabrazos Armrests	Abatible y extraíble. Folding and detachable.	
Reposapiés Footrests	Abatible interior-exterior. Altura regulable. Swing in and out. Height adjustable.	
Respaldo Backrest	No plegable. Bolsillo portaobjetos. Non-foldable. Back pocket.	
Accesorios Accessories	Anti-vuelcos. Reposapiés elevables. Porta oxígeno. Porta suero. Reposacabezas. Porta bastones. Anti-tipping wheels. Elevating footrests. Oxygen holder. Serum holder. Headrest. Cane holder.	